



Atomic Force Microscopy is a powerful technique for imaging, nanomanipulation, as platform for sensor work, nanolithography...

Conventional silicon or tungsten tips wear out quickly. CNT tip is robust, offers amazing resolution.



Simulated Mars dust




H. Dai




NASA Ames Research Center  
Ramsey Stevens, Lance Delzeit, Cattien Nguyen

## 280 nm Line/Space Array of Polymeric Resist on silicon Substrate




Silicon Tip

Multi-Walled Carbon  
Nanotube Tip




Regular Silicon Probe




MWNT Probe

**AFM images of a 280 nm line/space array of polymeric photoresist (300 nm thickness) on a silicon substrate acquired with a conventional silicon probe (left image) and a multiwalled carbon nanotube probe (right image).**


# AFM Imaging with Single Wall Nanotube Tips



2 nm thick Au on Mica



5 nm thick Ir on Mica



$\text{Si}_3\text{N}_4$  on Silicon substrate